

SOP for Determination of Bi-extensional Viscosity of Pounded Yam by Lubricated Squeezing Flow (LSF) Method

Biophysical Characterization of Quality Traits, WP2

Montpellier, France, January 2023


Santiago ARUFE VILAS, Centre de coopération Internationale en Recherche Agronomique pour le Développement (CIRAD), Montpellier, France

Oluwatoyin AYETIGBO, CIRAD, Montpellier, France

Romain DOMINGO, CIRAD, Montpellier, France

Christian MESTRES, CIRAD, Montpellier, France

Dominique DUFOUR, CIRAD, Montpellier, France (Validator)



This report has been written in the framework of RTBfoods project.

To be cited as:


Santiago ARUFE, Oluwatoyin AYETIGBO, Romain DOMINGO and Christian MESTRES (2023). *Standard Operating Protocol for Determination of Bi-extensional viscosity of Pounded yam by Lubricated Squeezing Flow (LSF) Method. Biophysical Characterization of Quality Traits, WP2*. Montpellier, France: RTBfoods Laboratory Standard Operating Procedure, 18 p. <https://doi.org/10.18167/agritrop/00686>

Ethics: The activities, which led to the production of this document, were assessed and approved by the CIRAD Ethics Committee (H2020 ethics self-assessment procedure). When relevant, samples were prepared according to good hygiene and manufacturing practices. When external participants were involved in an activity, they were priorly informed about the objective of the activity and explained that their participation was entirely voluntary, that they could stop the interview at any point and that their responses would be anonymous and securely stored by the research team for research purposes. Written consent (signature) was systematically sought from sensory panellists and from consumers participating in activities.

Acknowledgments: This work was supported by the RTBfoods project <https://rtbfoods.cirad.fr>, through a grant OPP1178942: Breeding RTB products for end user preferences (RTBfoods), to the French Agricultural Research Centre for International Development (CIRAD), Montpellier, France, by the Bill & Melinda Gates Foundation (BMGF). Authors thank the RTBfoods partner UAC-FSA and team for supplying the yam varieties and consumer tests data for the SOP from Benin.

Image cover page © LAJOUS P. for RTBfoods.



<div style="text-align: center;"> RTBfoods WP2: Biophysical Characterization of Quality Traits </div>		
SOP: Determination of Bi-extensional Viscosity of Pounded Yam by Lubricated Squeezing Flow (LSF) Method		
Date: 29/07/2022		Release: 1
<p>Written by: Santiago ARUFE VILAS</p> <p>For information on this SOP please contact:</p> <ul style="list-style-type: none"> • Oluwatoyin AYETIGBO, oluwatoyin.ayetigbo@cirad.fr • Santiago ARUFE VILAS, santiago.arufe_vilas@cirad.fr 		
<p style="text-align: center;">This document has been reviewed by:</p>		
Christian MESTRES (CIRAD) Oluwatoyin AYETIGBO (CIRAD) Romain DOMINGO (CIRAD) Dominique DUFOUR (CIRAD) Santiago ARUFE VILAS (CIRAD) Christian MESTRES (CIRAD) Oluwatoyin AYETIGBO (CIRAD) Santiago ARUFE VILAS (CIRAD)		03/01/2023 31/07/2022 01/08/2022 30/08/2022 21/12/2022 05/01/2023 05/01/2023 06/01/2023
<p style="text-align: center;">Final validation by:</p>		
Christian MESTRES (CIRAD)		06/01/2023

CONTENTS

Table of Contents

1	Scope and application.....	6
2	References.....	6
3	Principle	6
4	Reagents	7
5	Apparatus	7
6	Procedure	8
6.1	Calibration of texture analyser	8
6.2	Test conditions for measurement of BEV	8
6.3	Preparation of pounded yam dough and measurement.....	8
7	Expression of results.....	11
7.1	Obtained raw data	11
7.2	Method of calculation	11
7.2.1	Calculation.....	11
7.2.2	Expression of the results.....	12
7.3	Statistics	12
7.3.1	Analysis of Variance (ANOVA).....	12
8	Critical points or note on the procedure.....	12
9	Test report.....	13
9.1	ANOVA and repeatability of instrumental textural attribute.....	13
9.1.1	Bi-extensional viscosity (Pa·s) by Cultivar.....	13
9.1.2	Bi-extensional viscosity (Pa·s) by Cooking replicate	14
10	Appendices	15
10.1	Annex 1	15

ABSTRACT

This SOP describes the method for determining the bi-extensional viscosity (BEV) by application of lubricated squeezing flow (LSF) on pounded yam produced from different yam varieties using a texture analyser under fixed conditions. A full description of the sample preparation (pounded yam), conducting the experimental measurements, calculation using the obtained raw data and results expression are proposed. This SOP could be adapted for other pasty products such as Eba, Fufu or Matooke in order to obtain the corresponding bi-extensional viscosity that may be related to other sensory quality traits.

Key Words: bi-extensional viscosity, rheology, texture, compression, Hencky strain, strain rate.

1 SCOPE AND APPLICATION

This SOP describes the method for determining the bi-extensional viscosity (BEV) by application of lubricated squeezing flow (LSF) of pounded yam produced from varieties of yams using a texture analyser under fixed conditions.

2 REFERENCES

Otegbayo, B., Oroniran, O., Tanimola, A., Bolaji, O., Alamu, A., Mestres, C & Ayetigbo, O. (2022). Standard Operating Protocol for Textural Characterization of Pounded Yam. Biophysical Characterisation of Quality Traits, WP2. Iwo, Nigeria: RTBfoods Laboratory Standard Operating Procedure, 22 p. <https://doi.org/10.18167/agritrop/00613>

Ayetigbo, O., Domingo, R., Arufe, S., Mestres, C., & Akissoé, N., Otegbayo, B. (2022). Standard Operating Protocol for the Instrumental Determination of Extensibility of Pounded yam. Biophysical Characterization of Quality Traits, WP2. Montpellier, France: RTBfoods Laboratory Standard Operating Procedure, 18 p. <https://doi.org/10.18167/agritrop/00684>

Steffe, J. F. (1996). Rheological method in food process engineering. East Lansing, USA: Freeman Press, 2nd ed.

Van Vliet, T. (2008). Strain hardening as an indicator of bread-making performance: A review with discussion. Journal of Cereal Science 48, 1-9. DOI: <http://dx.doi.org/10.1016/j.jcs.2007.08.010>.

Launay G., Michon C. (2008) Biaxial extension of wheat flour dough: lubricated squeezing flow and stress relaxation properties. J Texture Stud 39:496–529. DOI: <https://doi.org/10.1111/j.1745-4603.2008.00156.x>

3 PRINCIPLE

The standard operating protocol (SOP) for determination of bi-extensional viscosity of pounded yam using lubricated squeezing flow measures the normal force applied to make the sample flow, the time and the distance of the texture analyse probe in order to calculate the Hencky strain and the strain rate during compression of pounded yam between two plates.

Squeezing flow between parallel plates can be achieved in many food rheology laboratories. When this deformation is executed between lubricated plates, biaxial extensional flow is achieved and bi-extensional viscosity can be calculated (Steffe,1996). In practice, three squeezing geometries may be used: (a) the upper and lower plates have the same diameter as the un-squeezed sample, (b) the upper compression plate and the sample have the same diameter, but the lower plate is much larger and (c) both the upper and lower compression plates are much larger than the sample. In the case of this SOP (a) geometry, that lead to a constant surface of contact between the sample and the plates was selected, Figure 1.




Figure 1. Lubricated squeezing flow with sample completely filling the gap between two parallel plates. R : radius of the plate, F_N : Normal force applied to the sample.

4 REAGENTS

Paraffin oil (CAS Number 8012-95-1)

5 APPARATUS

- Texture analyser (TA-XT Plus, Stable Micro Systems Ltd., Surrey, UK) with Exponent Software Interface.
- P/20P standard cylinder probe (20mm of diameter, Perspex, 10kg compression / 5kg tension)
- Yam Pounding machine (QZP/6000 model, Cheerfengly Ind. Co. Ltd, Taipei, Taiwan)
- Infra-red thermometer
- Chronometer
- Bowls
- Plastic moulds (14 mm height, 20 mm diameter, Teflon)
- Incubator (WTC binder, Tuttlingen, Germany)
- Sharp stainless-steel knives
- Cutting boards
- Mass balance ± 0.01 g (Precisa, Precisa Gravimetrics AG, Switzerland)
- Plastic spoon
- Plastic film roll

6 PROCEDURE

6.1 Calibration of texture analyser

The probe height and force are calibrated following the calibration instructions provided in the operation manual. The force calibration is done using a 2 kg standard weight. The test probe used is shown in Figure 2:



Figure 2. P/20P standard cylinder probe of 20mm diameter, Perspex (Stable microsystems).

6.2 Test conditions for measurement of BEV

Preliminary essays based on bibliography (Launay and Michon 2008, Van Vliet 2008) were carried out in order to find the best conditions for test measurement (please, see Annex 1 on the Appendix). After the evaluation of the preliminary data the selected test conditions were as follows:

- Test mode: compression
- Test probe: P/20P
- Load module: 5 kg
- Pre-test speed: 60 mm/min
- Test speed: 120 mm/min
- Post-test speed: 600 mm/min
- Trigger force: 0.25 N (25g)
- Target mode: distance
- Distance: 13 mm
- Sample dimensions: 20 mm of diameter and 14 mm of height

6.3 Preparation of pounded yam dough and measurement

The pounded yam dough is prepared following the RTBfoods SOP for pounded yam sample preparation recommended for texture analysis (Otegbayo et al., 2022) with some modifications as described below.

- Step 1: Two yam tubers are cut longitudinally into two halves and the proximal and distal sections are cut off from both ends (Figure 3).



Figure 3. Longitudinal halves of yam tubers after cutting proximal and distal sections.

- Step 2: From the two central sections, cubes of about 2 cm dimension are cut from both representative central sections (Figure 4).



Figure 4. Obtained yam cubes (≈ 2 cm side) from central section of the tuber after peeling and cutting.

- Step 3: 450 g of the yam cubes is weighed and placed in the cooking chamber of the yam pounding machine (Figure 5) above its receptacle containing 250 mL of water. The yam cubes are cooked for about 21 - 24 minutes (alarm sounds when cooking is finished), and pounded for 1.5 minutes in the pounding machine.



Figure 5. Yam Pounding machine (QZP/6000 model, Cheerfengly Ind. Co. Ltd, Taipei, Taiwan)

- Step 4: The yam dough is spooned out, 12 pre-lubricated moulds (20 mm of internal diameter, 14 mm of height) are filled; dough is cut out to fit their shape. The mould filled with dough sample is wrapped with plastic film (Figure 6).



Figure 6. Prepared samples.

- Step 5: These are then placed in an incubator at 55°C for 10 - 15 minutes to allow for stress relaxation and temperature equilibration.
- Step 6: The samples are demoulded, covered with paraffin oil and the compression test is applied with the aforementioned probe also covered with paraffin oil, for all replicates at constant temperature $\approx 30^{\circ}\text{C}$ (Figure 7).



Figure 7. Moulded samples lubricated with paraffin oil before and during measurement.

7 EXPRESSION OF RESULTS

7.1 Obtained raw data

After experimentation, raw data of normal force, F_N defined as the force applied perpendicularly to the sample, time (min) and distance (mm) are obtained during compression (total time is fixed to achieve compression until 1 mm).




Figure 8. Obtained raw data after measurement. Normal force (N) vs Distance (mm).

Once raw data are obtained, extensional viscosity can be calculated as described in Section 7.2.

7.2 Method of calculation

7.2.1 Calculation

Considering the configuration of the experiment, Figure 1, the **logarithmic strain** ε , also called *Hencky strain* of the sample is given by Eq. (1):

$$\varepsilon_H = -\ln \frac{h(t)}{h_0} \quad (1)$$

where h_0 is the initial sample height (14 mm) and $h(t)$ the height at time t . The biaxial (Hencky) strain \square_b of sample can be then calculated, Eq. (2):

$$\varepsilon_b = -\frac{1}{2} \ln \frac{h(t)}{h_0} \quad (2)$$

Once the biaxial (Hencky) strain \square_b is calculated, it can be related to the time in order to calculate the strain rate, $\dot{\varepsilon}_b$ (s^{-1}), Eq. (3).

$$\dot{\varepsilon}_b = \frac{d\varepsilon_b}{dt} \quad (3)$$

SOP: Determination of Bi-extensional Viscosity of Pounded Yam by Lubricated Squeezing Flow (LSF) Method

Date: 06/01/2023

Release: 1

The stress applied by the plates to compress the sample, σ (Pa) is calculated as the ratio between the normal force applied and the surface area of the plate:

$$\sigma = \frac{F_N}{S} \quad (4)$$

where F_N (N) corresponds to the normal force applied for dough deformation recorder by the software and S (m^2) is the surface of contact between the plates and the sample which is constant through the whole experiment ($S = \pi \cdot R^2 = \pi \cdot (10 \cdot 10^{-3})^2 = 3.14 \cdot 10^{-4} m^2$). Finally, $\dot{\epsilon}_b$ can be related to σ in order to obtain the bi-extensional viscosity, η_e (Pa·s), Eq. (5).

$$\eta_e = \frac{\sigma}{\dot{\epsilon}_b} \quad (5)$$

Once raw data are obtained, calculations are performed to obtain bi-extensional viscosity and can be facilitated in an Excel sheet (Microsoft Office Package) as follows:

Force [N]	Distance [mm]	Time [min]	h (t) [mm]	ϵ_H [-]	ϵ_B [-]	$\dot{\epsilon}_b$ [s^{-1}]	η_e [$N \cdot m^{-2}$]
				Eq. (1)	Eq. (2)	Eq. (3)	Eq. (4)

where blue cells correspond to raw data and grey cells to calculated data. Finally, a common η_e must be selected for all the studied samples. Once η_e value has been selected (*i.e.* $\eta_e = 1$), bi-extensional viscosity can be calculated using Eq. (5) and the corresponding values of σ and $\dot{\epsilon}_b$ associated to η_e .

7.2.2 Expression of the results

Results should be expressed using the International System of Units. The values of bi-extensional viscosity should be then expressed in Pa·s (or KPa·s). Moreover, the selected value of η_e (-) and the correspondent $\dot{\epsilon}_b$ must be clearly indicated (please, see Section 9).

7.3 Statistics

7.3.1 Analysis of Variance (ANOVA)

Experimental data were statistically analysed. The significant differences between means were identified by one-factor analysis of variance (ANOVA), followed by the Tukey test and considering significant P-values ≤ 0.05 using RStudio package (please, See Section 9).

8 CRITICAL POINTS OR NOTE ON THE PROCEDURE

- It must be ensured that the texture analyser is switched on at least 15 minutes prior to measurements and the texture analyser must be calibrated for force and distance on each experiment day prior to measurements.
- The dough must be relaxed for a period of 10-15 minutes at 50-55°C.
- Paraffin oil must be applied to each sample and to the probes between each measurement.
- Temperature of the dough must be monitored and measurement must be done at $\approx 30^\circ C$. Large variations in temperature may result in inaccurate results as the dough retrograde rapidly.

9 TEST REPORT

Table 1. Bi-extensional viscosity (KPa·s) of pounded yam cooking replicates obtained from different varieties at biaxial Hencky strain of 1, strain rate of 0.523 s^{-1} . (N: number of replicates, SD: standard deviation, SEM: standard error of the mean, CV: coefficient of variation of the mean).

Parameter	Cultivar	Cooking replicate	N	Mean	SD	SEM	CV
Bi-extensional viscosity (Pa.s)	Aga	1	2	106.6	52.5	37.1	49
		2	2	74.1	0.9	6.1	1
	Crachi	1	4	185.7	4.6	2.3	2
		2	4	166.7	10.0	4.9	6
	Dodo	1	8	50.7	6.9	2.4	14
		2	8	55.1	5.5	1.9	10
	Irindou	1	4	142.5	31.0	15.5	22
		2	4	116.9	52.7	26.4	45
	Laboko	1	4	286.0	17.2	8.6	6
		2	4	220.6	18.8	9.4	9
	Wete	1	4	139.9	13.8	6.9	10
		2	4	103.4	13.0	6.5	12

9.1 ANOVA and repeatability of instrumental textural attribute

9.1.1 Bi-extensional viscosity (Pa·s) by Cultivar




Figure 9. BEV of pounded yam varieties ($\varepsilon_b = 0.532 \text{ s}^{-1}$, $\varepsilon_B = 1$).

RTBfoods-WP2

SOP: Determination of Bi-extensional Viscosity of Pounded Yam by Lubricated Squeezing Flow (LSF) Method

Date: 06/01/2023

Release: 1

Table 2. One-way Analysis of Variance by cultivar.

	DF	Sum of Square	Mean Square	F value	Pr (>F)
Cultivar	5	$2.4 \cdot 10^{11}$	$4.8 \cdot 10^{10}$	68.54	$2.2 \cdot 10^{-16}$
Residuals	46	$3.2 \cdot 10^{10}$	$7.0 \cdot 10^8$		
C. Total	51	$2.7 \cdot 10^{11}$			

Table 3. Connecting Letters Report

	Aga	Crachi	Dodo	Irindou	Laboko	Wete
Level	AB	C	A	B	D	B

Table 4. Effect tests

	DF	Sum of Squares	Mean Square	F Ratio	Prob > F
Cultivar	5	$2.4 \cdot 10^{11}$	$4.8 \cdot 10^{10}$	107.4	$< 2 \cdot 10^{-16}$
Cooking replicate	1	$7.3 \cdot 10^9$	$7.3 \cdot 10^9$	16.3	$2.4 \cdot 10^{-4}$
Cultivar*Cooking replicate	5	$7.1 \cdot 10^9$	$1.4 \cdot 10^9$	3.2	$1.7 \cdot 10^{-2}$
Residuals	40	$1.8 \cdot 10^{10}$	$4.5 \cdot 10^8$		

9.1.2 Bi-extensional viscosity (Pa·s) by Cooking replicate




Figure 10. BEV of pounded yam varieties ($\dot{\epsilon}_b = 0.532 \text{ s}^{-1}$, $\epsilon_B = 1$) by cooking replicate.

Table 5. One-way Analysis of Variance by cooking replicate.

Source	DF	Sum of Squares	Mean Square	F Ratio	Prob > F
Cooking replicate	1	$7.3 \cdot 10^9$	$7.3 \cdot 10^9$	1.373	0.247
Residuals	50	$2.7 \cdot 10^{11}$	$5.3 \cdot 10^9$		
C. Total	51	$2.7 \cdot 10^{11}$			

RTBfoods-WP2	
SOP: Determination of Bi-extensional Viscosity of Pounded Yam by Lubricated Squeezing Flow (LSF) Method	
Date: 06/01/2023	Release: 1

Table 6. Connecting Letters Report

Cooking replicate	1	2
Level	A	A

The studied cultivars were significantly different in bi-extensional viscosity between one another, as shown in Table 3. The two-way ANOVA explained that both the effects of cultivar and cooking replicate on bi-extensional viscosity are significant, Table 4, as well as the significant effect of the interaction between cultivar and cooking replicate. The cultivar effect was the most significant on bi-extensional viscosity and the interaction effect is the least. The measuring bi-extensional viscosity of pounded yam by the proposed protocol may be considered as a repeatable, accurate and discriminatory protocol. Moreover, it was found a significant ($p \leq 0.10$) relationship between overall liking of the studied genotypes and BEV, Pearson $r=0.730$ (data not shown). This fact indicates that LSF method results are related to overall liking which makes this method a potential useful tool for screening of pounded yam.

10 APPENDICES

10.1 Annex 1

Preliminary essays with different compression speeds were carried out in order to determine the most discriminant compression speed in the case of pounded yam.

Four different genotypes, Aga (*D. Alata*), and Irindou, Laboko and Wete (*D. Rotundata*), were tested at 6, 12 and 120 mm/min compression speeds (corresponding strain rates of 0.026, 0.053, 0.532 s⁻¹, respectively) in order to obtain the corresponding bi-extensional viscosity, Figure 11.




Figure 11. Bi-extensional viscosity of pounded yam obtained from different genotypes at 0.026, 0.053 and 0.532 s⁻¹ strain rates.

Once bi-extensional viscosity of different genotypes at different strain rates was obtained, one-way ANOVA analysis was carried out in order to determine the best compression speed. Results, Table 7, shown that every tested compression speed lead to significant differences on the variance between the means of the populations, being the highest compression speed (associated to a strain rate of 0.532 s⁻¹) the one with the highest probability. Regarding the connecting letters reports after the application of Tukey test (Table 8), no difference of classification of groups for the studied population. Taking this into account, and in addition with the fact that the highest speed compression lead to the data obtaining in a shorter period of time, 120 mm/min was the compression speed selected for determination of bi-extensional viscosity of pounded yam.

Table 7. One-way Analysis of Variance by cultivar at a strain rate of 0.026 s⁻¹.

Strain rate (s ⁻¹)		DF	Sum of Square	Mean Square	F value	Pr (>F)
0.026	Cultivar	3	6.4·10 ¹²	2.1·10 ¹²	13.04	2.96·10 ⁻⁵
	Residuals	24	3.9·10 ¹²	1.6·10 ¹¹		
0.053	Cultivar	3	3.0·10 ¹²	9.9·10 ¹¹	10.75	1.11·10 ⁻⁴
	Residuals	24	2.2·10 ¹²	9.3·10 ¹⁰		
0.532	Cultivar	3	1.1·10 ¹¹	3.6·10 ¹⁰	28.25	4.81·10 ⁻⁸
	Residuals	24	3.1·10 ¹⁰	1.3·10 ⁹		

RTBfoods-WP2 SOP: Determination of Bi-extensional Viscosity of Pounded Yam by Lubricated Squeezing Flow (LSF) Method	
Date: 06/01/2023	Release: 1

Table 8. Connecting Letters Report as function of the used strain rate.

Strain rate (s ⁻¹)		Aga	Irindou	Laboko	Wete
Level	0.026	A	A	B	A
	0.053	A	A	B	A
	0.532	A	A	B	A

Results shown that the use a strain rate of 0.532 s⁻¹ (corresponding compression speed of 120mm/min) lead to the most discriminant results in terms of bi-extensional viscosity of the samples.



Institute: Cirad – UMR QualiSud
Address: C/O Cathy Méjean, TA-B95/15 - 73 rue Jean-François Breton - 34398 Montpellier Cedex 5 - France
Tel: +33 4 67 61 44 31
Email: rtbfoodspmu@cirad.fr
Website: <https://rtbfoods.cirad.fr/>