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C. CHAINE.

Description of an infrared oven for drying leaf samples

For leaf samples intended for chemical analysis to give a faithful picture of the nutrient status of the trees, a certain number of rules regarding the mode of sampling (« Conseils » No. 150) and of conditioning (« Conseils » No. 146) must be respected.

In the preparation of samples drying is of great importance. This operation must be done very soon after the samples are taken, and requires some form of desiccation apparatus (oven, heater).

If an adequate power supply is available, a practical method is to dry the samples in an electric oven heated by infrared rays, which is easy to construct. These « Conseils » describe an oven with four to six 250-watt infrared lamps and a capacity of 16 samples.

DESCRIPTION OF OVEN

The oven is made of aluminium sheet 2 mm thick; the sides are welded together. The front and back plates are two isosceles trapeziums, and are vertical (Fig. 1); each measures 123 cm at the broad base, 48 cm wide at the upper edge, and 53 cm high. The oven is 1 m long. The door in front opens on a horizontal hinge and is equipped with two hooks which fit over knobs on the side of the oven to keep the door open (Fig. 2). Its own weight holds the door closed, and it rests on angle-irons which serve as joints (Fig. 3).

This oven is designed for 6 screw-in, 250-watt infrared lamps (Fig. 4), or 1 500 watts in all. The power must be shared out by

distributing 500 watts (i. e. 2 lamps) per phase. The lampholders are immediately under the roof of the oven (Fig. 3). The distance between the heating elements and the trays of drying samples is 35 cm.

DESCRIPTION OF DRYING TRAYS

Each oven takes 4 drying trays measuring 23×92 cm, made of 1-mm aluminium sheet. Each tray is divided into 4 by 3 partitions (Fig. 5).

USE

A thermometer graduated 0-250 °C is placed horizontally between the drying trays; it is threaded through two pierced corks to isolate it from the hot metal. At the end of 1 hour's preheating the temperature should be stabilized at 105 °C maximum.

The oven takes 4 trays, say 16 leaf samples.

When placed in the drying trays, each sample should be accompanied by a Kraft paper label to identify it. Care must be taken to see that the labels do not get mixed up during handling. Drying is finished 2 hours after the samples are placed in the apparatus.

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