



Centre
de coopération
internationale
en recherche
agronomique
pour le
développement

Limits of accuracy and improvements on the HVI cotton fiber test

Cotton, from the plant to the final product

Jean-Paul Gurlot
Cirad-ca, Cotton Technology Laboratory


Plan of presentation

- **Introduction**
- **Different ways of fibers characterization**
- **A point about the standardization process**
- **How does work an HVI**
- **An example of relation between fiber and yarn quality**
- **Conclusions**

Cotton plant cycle




Cotton plant growth




From
Cirad
'Cotons'
software

Harvesting and ginning

Cotton bolls harvest




Cotton ginning




Saw

Statistics and economy




Cotton cycle* (1/3)




* From GOURLOT J.-P. et al 1999,
Recherche et développement en
technologie : mesurer et améliorer
la qualité des produits du cotonnier,
créer de nouveaux débouchés,
Agriculture et développement,
n°22, Juin 1999, ISSN 1249-9951,
pp. 90-113.

Commercial parameters
Length, Length unif, Strength
Grade, Micronaire
Average and variability


Cotton cycle (2/3)



Cotton cycle (3/3)




Fiber quality development



Fiber constitution

Constituants	Percentage
Cellulosis	95.0
Proteins	1.6
Waxes	0.9
Physiological sugars	0.3
Other	2.2


Fiber quality measurement



Fiber quality measurement

Could be measured

- Cellulose types
- Wax content
- Flexion / flexibility
- 'Curling' / crimp




Reversals **Spiral ...**

Plan of presentation

- Introduction
- **Different ways of fibers characterization**
- **A point about the standardization process**
- **How does work an HVI**
- **An example of relation between fiber and yarn quality**
- **Conclusions**

Fiber quality measurement

Why do we measure length ?



Define what type of final product will be made out of the fibers


Check the settings of ginning equipments

Define a commercial price of the fibers

Allow settings of the processing equipments

Fiber quality measurement

Manual :



pulling => commercial length

comb sorter => length diagrams (W, N)

=> ML, CV%, SFC Could be measured

Fiber quality measurement


Instrumental :

with classical instruments :

Fibrograph => SL2.5%, SL 50%, UR%

Almeter 101 => diagrams (W, N)

=> ML, CV%, SFC



Fiber quality measurement

Instrumental :

with classical instruments :


Fibrograph => SL2.5%, SL 50%, UR%

Almeter 101 => diagrams (W, N)

=> ML, CV%, SFC

with High Volume Instrument (HVI)

=> UHML, ML, UI%, SFI



Fiber quality measurement

Instrumental :

with classical instruments :

Fibrograph => SL2.5%, SL 50%, UR%

Almeter 101 => diagrams (W, N)

=> ML, CV%, SFC


with High Volume Instrument (HVI)

=> UHML, ML, UI%, SFI

with Advanced Fiber Information System (AFIS)


length (W, N)

=> ML, CV, UQL, SFC




Fiber quality measurement

Length histograms on 2 seeds (AFIS)





Fiber quality measurement

Length diagram (AFIS)




Fibre length

Part of the optical sensor



Fiber quality measurement


Length Fibrogram



$$\text{Uniformity Index \%} = \frac{\text{ML} \cdot 100}{\text{UHML}}$$

Fiber quality measurement

Length Fibrogram



$$UR\% = \frac{SL\ 50\ \% \cdot 100}{SL\ 2.5\ \%}$$



Fiber quality measurement

Length parameters

NOT FOR HVI


SL 2.5 ~ UHML ~ Pulling
22 – 40 mm

SL 50 % \leftrightarrow **ML**
10 – 18 mm 18 – 30 mm

UR % \leftrightarrow **UI %**
38 – 50 % 78 – 85 %

Fiber quality measurement

Why do we measure Fineness and Maturity ?



Predict the number of fibers per yarn cross-section


=> yarn strength and evenness

Chemical treatments and dye consumption

Dye uptake ability


Fiber quality measurement

Manual :



using microscope on longitudinal fibers

wall thickness using microscope on fibers sections




Fiber quality measurement

Manual :

with classical instruments

Fibronaire => Micronaire




Fiber quality measurement

Manual :

with classical instruments

Fibronaire => Micronaire



with Fineness Maturity Tester (FMT)

IM, MR, PM%, H, HS



Fiber quality measurement


Manual :

with classical instruments

Fibronaire => Micronaire

with Fineness Maturity Tester (FMT)

IM, MR, PM%, H, HS



with High Volume Instrument (HVI)

Micronaire



Fiber quality measurement

Manual :

with classical instruments


Fibronaire => Micronaire

with Fineness Maturity Tester (FMT)

IM, MR, PM%, H, HS

with High Volume Instrument (HVI)

Micronaire



wall thickness

with Advanced Fiber Information System (AFIS)

Distribution of Diameter,
Theta, MR, H



Fiber quality measurement

Recorded results:


IM : micronaire [2, 7]

MR : Maturity Ratio [0, 1.2]

PM% : Percent Mature fibers [0, 100]


H : Linear Fineness (mtex) [120, 350]

Hs : Standard Fineness (mtex) [120, 400]



$$H_s = \frac{H}{MR} \quad (\text{tex} = \text{grams}/1000 \text{ m})$$



Micronaire, maturity and fineness



=> Some combinaisons of MR and H correspond to close IM for very different fibres

Micronaire, maturity and fineness


IM = 4.1



MR = 1.04

H = 150

Hs = 144




MR = 0.67


H = 220

Hs = 328

Fiber quality measurement



Warp, IM = 6.0




Weft

IM = 2.4

IM = 3.5

IM = 4.7


Relation between IM, MR and Hs




$$MR^2 = (3.86 \times IM^2 + 18.16 \times IM + 13) / Hs$$

$$PM = (MR - 0.2) \times (1.565 - 0.471 \times MR) \times 100$$

Relation between IM, MR and Hs




Relation between IM, MR and Hs




Fiber quality measurement

Why do we measure color and trash ?



- To identify bales having homogeneous characteristics and group them per lot
- To avoid variations in color in raw and dyed fabrics
- To limit wastes during processing

Fiber quality measurement




Reflectance
Rd %

Yellowness
+ b


Image
analysis ...



Fiber quality measurement




HVI COLOR GRADES FOR AMERICAN UPLAND COTTON



Fiber quality measurement

Trashmeter image




Binarisation and threshold applied ...

Fiber quality measurement

Trashmeter image


8 pixels
considered
as trash
particle



Trash count, Trash area, Leaf, Leaf grade

Fiber quality measurement

Why do we measure strength ?



- To define what product can be made from these fibers
- To define a commercial price
- To predict yarn strength

Fiber quality measurement

Manual :

with classical instruments


Stelometer : T1 (cN/tex) and elongation E1 (%)

with High Volume Instrument (HVI)


Strength (cN/tex) and elongation Elong (%)

with devices measuring individual fibers


Breaking force and elongation



Fibrogram curve



Fibrogram curve after a break



Caractérisation de la fibre

Les principales caractéristiques de fibre peuvent être établies soit manuellement, soit sur appareils classiques ou soit sur chaînes HVI.

Seules les chaînes HVI :

permettent des mesures instrumentales, automatisées, rapides et intégrées des critères commerciaux,

et autorisent un classement balle à balle sur l'ensemble de ces critères

Fibre quality per spinning method

Rank	Ring spinning	Open end	Air-jet
1	Length	Strength	Length
2	Strength	Fineness	Trash
3	Fineness	Length	Fineness
4	Trash	Trash	Strength

Deussen, 1992

Fiber quality measurement



Example of ITC

(This display does not constitute any type of recommendation for this equipment, picture from an advertisement from Uster Technologies)

Example of HVI classification in Dumas (AK)




Fiber quality measurement

Actual major manufacturers of so-call HVI equipments (alphabetic order)

- Lintronic (Israel)
- Premier (India)
- Schaffner Technologies (USA)
- Zellweger Uster (USA)


Fiber quality measurement



From Hunter, 2000

Fiber quality measurement

High Volume Instrument (HVI) in the world




Source : HUNTER L., ITMF Committee, Brème, Mars 1994.

Plan of presentation

- Introduction
- Different ways of fibers characterization
- **A point about the standardization process**
- **How does work an HVI**
- **An example of relation between fiber and yarn quality**
- **Conclusions**

Normalization steps



Based on so-call 'reference methods'

Expected precision for international market

	Confidence interval
IM	(+/-) 0.1 unit
Length	(+/-) 0.02 inch (+/-) 0.51 mm
Length uniformity	(+/-) 1.5 %
Strength	(+/-) 1.5 cN/tex
Rd %	(+/-) 1 %
+b	(+/-) 0.5
Trash	(+/-) 0.1 %

From Sasser, 1992.



Reference cotton use in the world: 1991

Who	No HVI	Consumption			
		HVICC kg	Ratio kg/HVI	ICCS kg	Ratio kg/HVI
USDA	212	19214	91	3023	14.2
USA except USDA	91	828	9	49	0.5
Outside USA	318	192	0.6	1.8	0.004
Total	621				



Reference cotton use in the world: 1994

Who	No HVI	Consumption			
		HVICC kg	Ratio kg/HVI	ICCS kg	Ratio kg/HVI
USDA	x	14145 -21%	?	6803 +125%	?
USA except USDA	x	766 - 7%	?	127 +159%	?
Outside USA	x	245 +27%	?	42 +2200%	?
Total	883 +42%				




Existing reference cottons after 1998

ICCS	ICCS Mike only	HVICC
* Stelometer	HVI	HVI
* Carded cotton	Raw cotton	Raw cotton
* 3 types	6 types	2* 2 types
* Measurement of T1 E1 (1 type with SL%)	Mike	Length (UHM) UI% Strength




6 - International normalisation of the measurements

ICCS Mike only standards



Upper Half Mean Length Histogram
2356 Cottons from different origins




Length Maxi HVICC/Length Mini HVICC Tested cottons

6 - International normalisation of the measurements

Length of ICCS standards



Strength Histogram
2356 Cottons from different origins



St Maxi HVICC/St Mini HVICC Tested cottons



Existing reference cottons after 1998

UNIVERSAL
HIGH VOLUME
INSTRUMENT

Sept 2002
APPROVED

**LONG-STRONG
CALIBRATION
COTTON**

AMERICAN UPLAND


* 3 0 8 2 9 *

	STANDARD VALUE	STANDARD DEVIATION
Micronaire.....	3.99	0.06
Strength (g/tex).....	33.1	0.61
Upper Half Mean Length (in.).....	1.173	0.011
Uniformity Index (%).....	83.5	0.51


Micronaire Standard Deviation based on single specimen HVI testing. Standard Deviations of Strength, Length and Uniformity are based on two specimen HVI testing.

APPROXIMATE NET WT.
5 POUNDS

Existing reference cottons after 1998



HVI modules calibration : IM



About reference values establishment

Question : is there any drift in the establishment of Calibration Cotton reference data ?


- A set of calibration cotton = at least 2 cottons
- These can also be analyzed as samples

Calibration	Samples	Issue/Date
Set 1	Set 1, 2, 3, ...	1
Set 2	Set 1, 2, 3, ...	2
Set 3	Set 1, 2, 3, ...	3



Comparison between results => drift ?

Preliminary experiment



Drift and other troubles

- Observed Drift
- + Delay to access to recent generations of reference cottons


Consumption in 1991

=> Differential in reading levels



=> Lower accuracy in the yarn quality prediction from fiber quality measurement

Relations between yarn versus fiber quality



=> Use / consume reference materials !!!!




Economical incidence of quality on fiber exchange price



1 9 9 3




Evolution of premium/discounts for HVI strength over 5 years (P. SASSER EFS 1995)




Plan of presentation

- Introduction
- Different ways of fibers characterization
- A point about the standardization process
- **How does work an HVI**
- **An example of relation between fiber and yarn quality**
- **Conclusions**

Hardware design




Tenacity Material



**Force
sensor**

**Optical
sensor**


Tenacity Material




**Force
sensor**

**Traction
jaws**


Fibrogram curve




Fibrogram curve





Force / elongation curve




Fibrogram curve after a break




Fibrograms





Position of the break




Vue de face :



Vue de côté :




Vue de face :



pincettes avant

Vue de côté :



Rheology

Constraint


**Elasticity
of
reorganization**

Elasticity

Slippage

Plasticity


Deformation




Tenacity or strength

Results and conclusions

HVI




Single fibers




Tenacity $\cong F / S$ should be constant

Optical sensor



Possible effects of fiber properties (example on strength)



Attention: the mass of broken fibers is estimated from this signal for the calculation of strength

- 1-Sensor linearity
- 2-Shape factor
 - micronaire
 - scouring / swelling
- 3-Maturity and fineness distributions
- 4-Color
- 5-Length distribution
- 6-Ambiant conditions

Ambiant conditions

Should be at any time:

- 21 °C +/- 1 °C
- 65 % Relative Humidity (RH%) +/- 2 %



Ambient conditions


Material and method

Optical sensor

Brush

Hypothesis : ambient conditions should be stable around the sensor

- Use of temperature probes
- Use of Comb




Ejection

Fibrosampler

Air ejection OFF
ON



Ambient conditions

Results and conclusions : air suction OFF



Ambiant conditions


Results and conclusions : air suction ON



Moisture content (% w/w)

Source: USDA

⇒ More than 1 cN/tex
difference !
(Sasser 1990)






Strength stability : HVICC bale 27985

Measuring unit covered vs uncovered

Motion Control 3500







Strength stability : HVICC bale 28484

Measuring unit covered vs uncovered

Motion Control 3500




Other possible bias in HVI measurement



Mass of this deposit = 0.0175 g

Sources of variability in the results



Variety
Growing conditions
(fertilizer, insects)



Plant to plant
Picking technique
Farm size



Seed cotton preparation
Ginning technique (R/S)
Lint cleaning




Number of samples / bale
Number of bales / lot
Method of sampling

**Precision
Accuracy
Repeatability
Reproducibility**



RH conditions,
HVI calibration
Nb tests / sample

Confidence Intervals (research samples in specific sampling conditions)




Upland Saw
ginned
samples

Barbadosse
cottons

Plan of presentation

- Introduction
- Different ways of fibers characterization
- A point about the standardization process
- How does work an HVI
- **An example of relation between fiber and yarn quality**
- **Conclusions**

Relations between yarn versus fiber quality




Relations between yarn versus fiber quality

	Coton A	Coton B
ML (mm)	22.6	24.2
UHML (mm)	28.4	28.9
UI (%)	79.6	83.7
Strength (cN/tex)	25	30.3
Elong (%)	5.0	5.7
IM	2.9	3.8
MR	0.67	0.90
PM (%)	58.5	79.6
H (mtex)	143	155
HS (mtex)	214	173
Rd (%)	69.7	72.8
+b	11.9	11.3


Relations between yarn versus fiber quality

CRL yarn tenacity tests




Relations between yarn versus fiber quality

CRL yarn tenacity tests




Relations between yarn versus fiber quality

CRL yarn tenacity tests (RS)




Relations between yarn versus fiber quality

Evenness tester UT3 : thick places (RS)




Relations between yarn versus fiber quality

Evenness tester UT3 : thin places (RS)




Relations between yarn versus fiber quality

Eveness tester UT3 : neps places (RS)




Relations between yarn versus fiber quality

Eveness tester UT3 : hairiness (RS)




Relations between yarn versus fiber quality

Evenness tester UT3 : CV% (RS)



Relations between yarn versus fiber quality




Effect of fibre parameters on yarn resistance

- 191 cottons from various origins
- Fibre analysis + spinning RS 20tex


$$\text{Ten Fil} = 0.44 \text{ TenHVI} - 0.0016 \text{ H} + 2.58 \text{ MR} + 0.33 \text{ UI} - 27.03$$
$$R^2 = 0.76 ***$$

Fibre strength and UI vs yarn strength RS 20 tex


UI%




Fibre strength and UI vs yarn strength RS 20 tex




Fibre strength and MR vs yarn strenght RS 20 tex




Fibre strength and MR vs yarn strenght RS 20 tex




Fibre strength and H vs yarn strenght RS 20 tex



Fibre strength and H vs yarn strenght RS 20 tex



Fibres characteristics vs yarn strenght RS 20 tex



Fibres characteristics vs yarn evenness RS 20 tex

- 30 cottons
- Fibres characterization
- Spinning OE 20, 27 and 37 tex
- Spinning RS 20, 27 and 37 tex

Correlations coefficients between fibres characteristics and OE yarn evenness

		ML	UHML	UI	ST	EL
Thin pl. (FIN)	20 tex	-0.55	-0.58	-0.33	-0.62	-0.43
Thick pl. (GRO)	“	-0.34	-0.37	-0.12	-0.34	-0.20
Neps (NEP)	“	-0.25	-0.28	-0.10	-0.26	-0.22
CV evenness (CVR)	“	-0.62	-0.63	-0.48	-0.60	-0.49
Hairiness (PIL)	“	-0.63	-0.61	-0.60	-0.55	-0.59
FIN	27 tex	-0.54	-0.57	-0.33	-0.52	-0.41
GRO	“	-0.32	-0.35	-0.12	-0.29	-0.25
NEP	“	-0.29	-0.32	-0.10	-0.28	-0.27
CVR	“	-0.39	-0.40	-0.27	-0.27	-0.34
PIL	“	-0.57	-0.56	-0.52	-0.52	-0.41
FIN	37 tex	-0.43	-0.44	-0.32	-0.34	-0.30
GRO	“	-0.39	-0.41	-0.23	-0.28	-0.34
NEP	“	-0.25	-0.29	-0.03	-0.21	-0.17
CVR	“	-0.29	-0.28	-0.29	-0.21	-0.26
PIL	“	-0.65	-0.64	-0.59	-0.49	-0.62

Correlations coefficients between fibres characteristics and RS yarn evenness

		ML	UHML	UI	ST	EL
Thin pl. (FIN)	20 tex	-0.75	-0.70	-0.86	-0.68	-0.71
Thick pl. (GRO)	“	-0.85	-0.81	-0.87	-0.78	-0.81
Neps (NEP)	“	-0.60	-0.61	-0.42	-0.51	-0.60
CV evenness (CVR)	“	-0.89	-0.86	-0.89	-0.82	-0.83
Hairiness (PIL)	“	-0.72	-0.66	-0.86	-0.64	-0.73
FIN	27 tex	-0.71	-0.66	-0.85	-0.64	-0.67
GRO	“	-0.82	-0.80	-0.79	-0.75	-0.79
NEP	“	-0.49	-0.53	-0.24	-0.41	-0.48
CVR	“	-0.90	-0.87	-0.90	-0.84	-0.85
PIL	“	-0.68	-0.62	-0.84	-0.64	-0.72
FIN	37 tex	-0.66	-0.60	-0.82	-0.58	-0.60
GRO	“	-0.81	-0.79	-0.74	-0.73	-0.75
NEP	“	-0.47	-0.51	-0.20	-0.40	-0.43
CVR	“	-0.90	-0.86	-0.90	-0.83	-0.84
PIL	“	-0.71	-0.66	-0.86	-0.64	-0.75

Plan of presentation

- Introduction
- Different ways of fibers characterization
- A point about the standardization process
- How does work an HVI
- An example of relation between fiber and yarn quality
- **Conclusions**

Conclusion

HVI cotton fiber measurements may be used to :

- **Commercially characterize cotton fibers properties**
 - standardization ongoing for all measured parameters
 - future evolutions to integrate new properties characterization
 - that may induce new rules in the trade

Conclusion

HVI cotton fiber measurements may be used to :

- **Arrange laydowns to stabilize or control :**
 - mean values
 - variability around those mean values according to production means (from field to ginning mill), sampling procedures (from ginning to spinning mills) ...

Conclusion

HVI cotton fiber measurements may be used to :

- **Predict the fiber behavior in the processing steps both in terms of :**
 - quality
 - productivity

Conclusion

HVI cotton fiber measurements may be used to :

- **Control,**
- **Check,**
- **And set spinning machineries**

to get the highest yarn quality as demanded by the market

Conclusion

HVI cotton fiber measurements may be used to :

- **Breed new varieties depending on the improvements made in the transformation stages.**

depending on commercially recognised
characterization

What Cirad does and recommends

- 1) Apparatus calibration to insure a proper reading level.
- 2) Check and set-up of procedures to warrant proper precision and accuracy levels.
- 3) Check of the results through a participation to periodical international round tests.
- 4) Check the precision in classing routine.

TO GET

- Homogeneous results on the cotton market
- Limited number of claims.



Centre
de coopération
internationale
en recherche
agronomique
pour le
développement

Thank you
for your attention

