



Root system structure

A. Audebert

2nd Global Rice Phenotyping Network Workshop
IRRI 22-24 November 2012



Rhizoscope :

A phenotyping system for root dynamics



Objective : Adult plants phenotyping



Principle : Mixed hydroponic system with glass beads as inert substrate

- Transparent Plexiglas boxes (50 cm x 20 cm x 2 cm) with a grid of nails
- Filled with 1.5 mm diameter glass beads = inert substrate providing mechanical resistance.
- Automatic control of pH and temperature of the nutrient solution
- Medium throughput system: 4 tanks holding each 48 plants (192 plants)
- Plants up to 30 days old





What traits can be monitored?

- Manually or through image analysis:
 - Shoot traits
 - Leaf length; number of tillers per plant; leaf number; shoot biomass;
 - Root traits:
 - Maximum crown root depth; root mass at different depths; crown root number at depth; root cone angle
 - Root /shoot ratio; root mass per tiller
- Fluorescence microscope on crown root tips
 - Root diameter; stele area; nb of xylem vessels; diameter of xylem vessels; nb of mesodermis layers




Phenotyping experiments

- 2011 and 2012
- Orytage and Grips funds
- Japonica and Indica panels (172 varieties phenotyped)
- 2 checks varieties (IR64 and Azucena) replicated in each tank
- 2 replications in time
- 30 days growth
- Traits
 - Shoot biomass
 - Root cone angle
 - $\alpha = 180 - (\beta_1 + \beta_2)$
 - Maximum root depth
 - Root nb at 30 cm depth
 - Root biomass
 - 0-20
 - 20-30
 - 30-50
 - Anatomy
 - Root tip




Rhizoscope platform : Some results

(indica panel, 200 acc.)






Rhizoscope platform : Some results

(indica panel, 200 acc.)






Rhizoscope platform : Some results

Indica panel



japonica panel




Sub-population 1 : Indicas from Madagascar

How explain the japonicas diversity ?

Available information for theses varieties: Type (traditional or improved);
Geographical origin (large region); Genetic group


Diversity tree
established from 25
Markers



Sub-population 2 : Japonicas from Africa and South America

Sub-population 1 : Japonicas from Asia


PCA analysis



Next step

- Association study
 - With the SNPs genotyping data
- Improved the picture analysis
 - System
 - Software
- Growth condition
 - Nutrient
 - Salt, N, P
 - Abiotic
 - pH, temp, O₂, impedance
- Improved beads capillarity





Thank you !