

20-21-22 MARS  
PALAIS DES CONGRÈS  
PORTE MAILLOT - PARIS  
2013

IRC

INTERNATIONAL  
RUBBER  
CONFERENCE

**IRC**  
International Rubber  
Conference Organisation

AFICEP

dgcis

direction générale de la compétitivité  
de l'industrie et des services



**Programme**

[www.irc2013.com](http://www.irc2013.com)

## AFTERNOON : PARALLEL SESSIONS / APRÈS-MIDI : SESSIONS PARALLÈLES

	THEME 2 PROCESSING, OPTIMIZATION OF PROCESSES TECHNIQUES DE TRANSFORMATION ET OPTIMISATION DES PROCÉDES	THEME 1 CHEMICALS ASPECTS AND RAW MATERIALS ASPECTS CHIMIQUES, MATIÈRES PREMIÈRES ET FORMULATION
14.00 - 14.40	<b>Guest speaker / Conférencier invité</b> Prof.Dr.-Ing.Andreas LIMPER, Member of the Management Board, HARBURG-FREUDENBERGER Maschinenbau GmbH, (Germany) New angles in mixing of rubber compounds	<b>Guest speaker / Conférencier invité</b> Dr. Luca GIANNINI Nanofiller Project Leader, PIRELLI TYRE, (Italia) A Rubber Vision of High Aspect Ratio Nanofillers
14.40- 15.00	Rheological behavior of NR by LAOS analysis in relation with postharvest and processing parameters <b>F. DEME - LRCCP</b> <b>J. SAINTE-BEUVE - CIRAD (France)</b>	New material solutions for harsher service conditions of automotive hose <b>S. BOUVIER - DUPONT PERFORMANCE POLYMERS</b> (Suisse)
15.00- 15.20	Efficience énergétique & réduction des pertes matière lors de la fabrication des articles moulés en caoutchouc <b>J-L. MAIRE - DESMA ELASTOMERTECHNIK (Germany)</b>	Synthesis of telechelic polyisoprenic oligomers by cross- metathesis reaction from natural or synthetic rubber <b>A. NOURRY - UNIVERSITE DU MAINE (France)</b>
15.20 - 15.40	Industrial compoundings of elastomers with nanofillers <b>M. CHARMAN - EMAC (France)</b>	Modified SSBR for Silica & Carbon Black Containing Tires <b>S.THIELE - STYRON EUROPE GMBH (Germany)</b>
15.40 - 16.00	Continuous rubber recycling using a co-rotating twin screw extruder <b>A. GALLO - F.LLI MARIS (Italia)</b>	Noxite® High Performance High Temperature Acrylate Rubbers (HT-ACM) Improved media resistance in fuels and new highly addivated engine and transmission oils - A comprehensive study <b>K. ZOUMIS - UNIMATEC CHEMICALS EUROPE</b> (Germany)
16.00 - 16.20	Pause and visit of the exhibition	Novel polyurethane foams based on hydroxytelechelic oligomers of natural rubber <b>T. K. N. TRAN - UNIVERSITE DU MAINE - (France)</b>
16.20 - 16.40	Elastomer network investigation by solid NMR <b>B. GABRIELLE - E. GOMEZ - HUTCHINSON (France)</b>	Pause and visit of the exhibition
16.40 - 17.00	Rubber to metal bonded assemblies Characterization of the adhesion for different formulas of elastomers, metals and bonding agents. Correlation between kinetic crosslinking of rubber compounds and bonding agents <b>F. BRUNO - LRCCP (France)</b>	Identification des facteurs de variabilité du caoutchouc naturel pour développer des nouveaux grades à variabilité réduite et contrôlée destinés aux applications antivibratoires - Projet collaboratif CANAOPT <b>J. SAINTE-BEUVE - CIRAD</b> <b>J-F. PILARD - UNIVERSITE DU MAINE (France)</b>
17.00 - 17.20	How to reduce fatigue life dispersion starting from the study of 3D distribution of ZnO agglomerates by computed X-ray microtomography <b>K. LE GORJU - HUTCHINSON (France)</b>	HT-ACM XP - The Winning Combination of Processing and Performance for High-Temperature, Oil Resistant Hose <b>P. ABRAHAM - ZEON CHEMICALS EUROPE</b> (England)
17.20 - 17.40	Electron Beam Irradiated Natural Rubber-Polyaniline Dodecylbenzenesulfonate Blends <b>Y. KOK-CHONG - MALAYSIAN RUBBER BOARD</b> (Malaysia)	GUAYULE / HEVEA latex glove comparison <b>S. PALU - CIRAD (France)</b> <b>M.DORGET-CTTM (France)</b> <b>D.GUERIN-PIERCAN (France)</b>

**ansm**

Agence nationale de sécurité du médicament  
et des produits de santé



**cirad**

LA RECHERCHE AGRONOMIQUE  
POUR LE DÉVELOPPEMENT



**CTTM**  
CENTRE DE TRANSFERT  
DE TECHNOLOGIE DU MANS

# **GUAYULE / HEVEA**

## **latex glove comparison**

**M. DORGET,  
A.AMOR,  
S. PALU,  
D.PIOCH  
E.TARDAN,  
D. GUERIN,  
C.MOURTON-GILLES,**

**CTTM  
CTTM/CIRAD  
CIRAD  
CIRAD  
CIRAD  
PIERCAN  
ANSM**

[www.cttm-lemans.com](http://www.cttm-lemans.com)





**Hevea plantation,**



**Guayule field , Spain (El Molinar)**

## WHY AN ALTERNATIVE SOURCE OF NR SUPPLY ?

- **HEVEA**, only commercial source of NR ( 93% of world production in Asia),
- Growing demand from emerging countries (China, India, Brazil),
- NR prices have rocketed upward,
- Threat by *Microcyclus ulei* (SALB) (only in South America). The risk to spread in Asia/Africa exists. When ? How ?
- Price of NR and SR linked with oil (Brent), volatile price



## **WHY AN ALTERNATIVE SOURCE OF NR SUPPLY ?**

- **Proteins in Hevea latex cause life-threatening, IgE-latex allergy. Guayule latex hypoallergenic,**
- **Tendency for replacement rubber plantations by palm oil plantations,**
- **World climatic changes, more frequent weather woes.**
- **Rubber tapping laborious, social aspects (NGO), Guayule can be mechanized,**
- **Hevea, Guayule, TKS, Polyisoprene cis 1-4, High Mw,**



# LATEX PROPERTIES

	HEVEA latex	Commercial GUAYULE latex
Solid content (%)	<b>61.4</b>	<b>55.6</b>
Viscosity (Cp)	<b>48</b>	<b>53</b>
pH	<b>9.6</b>	<b>10.9</b>



**Similar wet characteristics**



# DRY CHARACTERISATION

drying procedure :


- 2 weeks at ambient condition
- 2h. at 110°C



1. concentration
2. deformation
3. coalescence




honey comb structure




No additive, no vulcanisation



# DRY CHARACTERISATION-STRAIN



# DRY CHARACTERISATION- LINEAR DOMAIN



- HEVEA film behaves like a “solid” film as far as 200°C
- GUAYULE film behaves like a “liquid” film since 100°C


 **Very different dry characteristics**

# GLOVES PRODUCTION


1. coagulant surface treatment
2. latex surface coagulation
3. compaction



free particles



Coagulant surface treated



coagulated particles



surface coagulation process

with additives and vulcanisation

# GLOVES PRODUCTION



quality progress after formulation and process adaptations






# GLOVES PROPERTIES

	HEVEA LATEX	COMMERCIAL GUAYULE LATEX	EU-PEARLS GUAYULE LP LATEX	EU-PEARLS GUAYULE HP LATEX
Solid content (%)	<b>61.4</b>	<b>55.6</b>	<b>37.7</b>	<b>15.5</b>
Sleeve thickness (mm)	<b>0.32</b>	<b>0.28</b>	<b>0.14</b>	-
Hand thickness (mm)	<b>0.37</b>	<b>0.31</b>	<b>0.14</b>	-

➔ **GUAYULE  $\approx$  HEVEA**  
gloves thickness if solid  
contents (DRC)  
are the same



# GLOVES PROPERTIES

	HEVEA LATEX	COMMERCIAL GUAYULE LATEX	EU-PEARLS GUAYULE LP LATEX	EU-PEARLS GUAYULE HP LATEX
<b>Vulcanising dispersion ratio</b>	<b>27</b>	<b>54</b>	<b>54</b>	<b>54</b>
<b>Stress at break (Mpa)</b>	<b>17</b>	<b>7.0</b>	<b>13.3</b>	<b>-</b>
<b>Strain at break (%)</b>	<b>810</b>	<b>860</b>	<b>808</b>	<b>-</b>

 **GUAYULE and HEVEA gloves mechanical properties are similar after slight formulation and process adaptations**

# Analysis of protein extracts from films and gloves produced with hevea and guayule latex

**C. MOURTON-GILLES, M. DANIAU, G. OLIVIER  
A. SOUVERAIN, S. PALU, D. PIOCH, E. TARDAN  
M. ROSSIGNOL, S. HEM**

**ANSM  
CIRAD  
INRA-Protéome**



## Scientific context

**Latex allergy is due to the présence of hevea proteins** in rubber goods (cleaning gloves, nipples, toys,...) and medical devices (gloves, catheter, syringe,...).

### **Latex allergy is a real health problem**

Up to 17% of health care workers suffer from type I hevea allergy

Up to 16% of intra-operative shock are due to this allergy

Presence of cross allergy with some fruits (banana, kiwi,...)

A solution to this problem could be to develop devices with a non allergenic latex : **GUAYULE LATEX**

## Objective

Collaborative study between CIRAD and National Health product Agency

To compare protein content of Guayule and Hevea device extracts






## DETERMINATION PROTEINS CONTENT ON FILMS AND GLOVES EXTRACT

**-Extraction of proteins**  
*AFNOR EN 453-3 Modified*


**-Dosage proteins**  
*Lowry Kit Biorad DC*



- Proteins quantity in guayule films are 115 times lower than in hevea films and 9 times than in medical gloves (1)
- Gloves made with commercial guayule latex contain less proteins than hevea prototype gloves (2)

## DETERMINATION OF ALLERGENS CONTENT FILMS AND GLOVES

-Dosage of allergens  
*Elisa sandwich*  
*FITKit ICOSAGEN*  
(Acm, ABTS)




- Allergen Hevb1 is present in the hevea films but in small amount in guayule extracts from films and gloves

*Commercial guayule latex was provided by YULEX for the EU- PEARLS project*

## DETERMINATION OF ALLERGENS CONTENT FILMS AND GLOVES

-Dosage of allergens  
*Elisa sandwich*  
*FITKit ICOSAGEN*  
(Acm, ABTS)




- Allergen Hevb 3 is present in the hevea films and in small amount in medical and guayule gloves made with commercial latex

*Commercial guayule latex was provided by YULEX for the EU- PEARLS project*

## DETERMINATION OF ALLERGENS CONTENT FILMS AND GLOVES

-Dosage of allergens  
*Elisa sandwich*  
*FITKit ICOSAGEN*  
(Acm, ABTS)




➤ Allergen Hevb 5 is present in the hevea films

Commercial guayule latex was provided by YULEX for the EU- PEARLS project



## DETERMINATION OF ALLERGENS CONTENT FILMS AND GLOVES

-Dosage of allergens  
*Elisa sandwich*  
*FitKit ICOSAGEN*  
(Acm, ABTS)



- Allergen Hevb 6.01 is present in the hevea films but also in guayule gloves film in smaller quantity

*Commercial guayule latex was provided by YULEX for the EU- PEARLS project*

## Identification of proteins in guayule films extracts

**Using Western Block and labelling with anti-hevea Mabs (Hevb1,Hevb3,Hevb5 and Hevb6.02) and using LC ESI MS/MS mass spectrometry with INRA:**

➤ **Some hevea allergens were identified in guayule films**

- Rubber elongation factor Hevb1
- Glucan endo1,3 beta glucosidase: Hev b2
- Hevamin A : Hev b14
- Pro-Hevein: Hev b6.01

- **Proteins in guayule films extracts, Hevb1 and Hevb6.01 were detected.**
- **Allergen Hevb 6.01 is present in hevea films and gloves but also in guayule film and gloves**
- **Protein quantity in guayule films is about 100 times lower than in hevea films**
- **Only Hevb1 and Hevb6.01 were detected in guayule films using Fit Kits**

## CONCLUSIONS

- **For wet characterization, HEVEA and GUAYULE latex are very similar,**
- **For dry characterization, HEVEA latex behaves as “solid” film , GUAYULE latex as “liquid” film. With higher molar mass ( $M_w$ ) of the polyisoprene, GUAYULE films have similar or even better behavior than HEVEA films.**
- **For gloves production, slight formulation and process adaptations bring close mechanical properties and process behaviors for both latex.**
- **Preliminary results on latex allergenic of Hevea & Guayule latex, show that proteins concentration of guayule film is 100 fold lower than for HL film.**
- **Guayule latex has a lower content of allergenic proteins than Hevea latex.**
- **Guayule latex is an alternative source very useful in medical applications to avoid allergenic reactions of patients for Type I latex allergy**



**THANK YOU**

