



New alternative spectroscopic method for the detection of foodborne pathogens

Samira SARTER , Philippe DANIEL

CIRAD -UMR Qualisud

Institut des Molécules et des Matériaux du Mans

IMMM UMR CNRS 6283

Food safety risks

Food safety risks

Salmonella spp.

- Raw meat sold in market: Porc 39-64%; chicken 42-49-53%; beef 62%
- Resistance in meat: Porc 50-73% ; Chicken 45%
 - Tetracycline, sulphonamide, streptomycin, ampicillin, chloramphenicol, trimethoprim, nalidixic acid
- Multiresistance : 21-56% of isolates
 - 7-9 antibiotics: 15% / 10-13 antibiotics: 8%
- Multiresistant *Salmonella* from food or food-producing animals are common in different countries:
 - Malaysia 49-75% (n=88)
 - Thailand 44-66% (n=342)
 - Vietnam 21-56% (n=180)



Food safety risks

Listeria monocytogenes

- EU rejections: Filet Pangasius (8 notifications 2010; 17 en 2009)

Campylobacter spp.

- Chicken sold in market: 15.3%
- Chicken : 95% of strains are resistant to fluoroquinolones (critical AB)

Escherichia coli : a reservoir

- Resistance: 84% of isolates of beef, poultry, porc
- Resistance to fluoroquinolones: 16-21% of isolates, mainly in chicken samples (52-63%)
- Multiresistant *E. coli* (n=99) in raw meat:
 - 89.5% in chicken meat
 - 95% in chicken faeces
 - 75% in pork meat isolates

Food safety risks

- **Food Safety Objectives:** "the maximum frequency and/or concentration of a hazard in a food at the time of consumption that provides or contributes to the appropriate level of protection (ALOP)".
- To ensure that an FSO is met, it is required to set **Performance Objectives** which correspond to the levels that must be met at earlier steps in the food chain before consumption.
- FSOs and POs must be achievable by the application of good practices (GAP, GHP, GMP) and HACCP
- **Microbiological Criteria** can be used to define the microbiological quality of raw materials, food ingredients, and end-products at any stage in the food chain.

**Need for accurate, rapid and sensitive methods for detection
and quantification of microbial hazards**

Standard methods for pathogen identification

AFNOR ISO 6579:2002

Identification of *Salmonella* spp

Incubation

2 - 4 days

Identification

Time depending on method

Many hours

Phenotypic methods

**Immunological methods
(ELISA)**

**Molecular methods
(PCR)**

Biochemical methods

25g of sample




**Pre enrichment
Incubation in BPW**




**Selective enrichment
RVS + MKTTn**



**Isolement
XLD + XLT4**



**Incubation
Agar plate**



Applications of Raman spectroscopy to bacteria


Principles of Raman spectroscopy

monochromatic visible radiation : Laser ω_0, λ_0

Interaction with a sample




Raman effect gives the vibrational signature of any kind of materials



Scattered radiations

Inelastic process




Sir Chandrasekhara
Venkata RAMAN
1888-1970

Advantages of the technics:

- Fingerprint technics
- No preparation of the sample
- Non invasive technics
- Non destructive technics
- Qualitative or quantitative

Number of publications related to Raman scattering and bacteria


Source : ISI Web of Science – January 2014 –
Key words: Raman, bacter*




Raman study of bacteria



- Investigation of microcolonies and characterisation of heterogeneity



- Single-cell analysis of bacteria



Raman study of bacteria

Interprétation of the spectrum: ***fingerprint technique***


Exemple of E-coli

Lipids

Carbohydrates

Proteins

Nucleic acids



507 : Carbohydrate C-O-C

652 : Tyrosine (Acide Aminé)

727 : Adénine (ADN)

872 : Tyrosine (Acide Aminé)

955: Lipides

1037 : Lipides

1240 : amide III

1323 : $\delta(\text{CH}_2)$

1377 : Symm Stretch (CON-), $\delta(\text{CH}_2)$


1464 : mono-oligosaccharides

1580 : ADN


1771 : Ester

Raman study of bacteria

Allow to distinguish between types of bacteria




0,5 to 3 μm




Type de liaison	Raman shift cm^{-1}
Carbohydrates	460 ; 590
Acides nucléiques	770
Tyrosine	850
Phénylalanine	980 ; 1002
Amide III	1100
Lipides	1240
Amide II	1440
Amide I	1630 ; 1705
Lipides insaturés	1630 ; 1705

Gram +	Gram-
<i>Bacillus subtilis</i>	<i>Escherichia coli</i>
<i>Staphylococcus</i>	<i>Salmonella</i>
<i>Streptococcus</i>	<i>Pseudomonas</i>





Ward's algorithm
Gammes spectrales 400-1800 cm^{-1}

Raman study of bacteria by Raman spectroscopy vs growth phases



— Latence phase
— Exponential phase
— Stationary phase



L. Bendriaa, PhD Thesis , 2005

« Rather easy » distinction between young bacteria and old bacteria


Frequency range used for classification:
1450-1750 cm⁻¹

Functionalized surfaces for detection of pathogenic microorganisms

Alternative method

Biosensor based on a « **double check procedure** » :

- (1) Specific capture of microorganisms
- (2) Recognition by Raman spectroscopy




Specific functionalized
surface




Raman spectroscopy
analysis

Identification via
spectra recognition




Statistical data
analysis




A result of presence/ absence of pathogens in less than 24h


Exemple: Gold surface functionalisation with parabenzenesulfonyl chloride







Synthesis of specific surfaces of gold with chemical modifications



Protein A




Antibody




Antibody – antigen specific recognition


Capture of the microorganism




Protein A




Antibody



Raman characterization




Fluorescence image




Capture of the microorganism

➔ Evidence of the last step of the process




Anti-IgG
(1,07g/l)




➤ Functionalization procedure also possible on other type of substrate :

- Polyethylene treated by plasma
- Functionalized Polyurethane
- Systems including nanoparticles (magnetic, silver, gold: SERS effect)

QCM monitoring




Raman characterization



Raman spectra (785 nm, 10 mW) of
Salmonella immobilized on
functionalised Au surface

Conclusion and Perspectives

- Develop a detection kit based on Raman spectroscopy for specific pathogens in food (model and food matrix)
- Target specific resistant bacteria, and try to explore the mechanisms of actions (critical antibiotics)
- Screening of resistant strains along the food chain/environment
- Research at the interface between physics and chemistry of materials



Institute for Molecules and Materials of Le Mans

Department of solid state physics:

- **Physics of advanced materials, Nanomaterials, Surface functionalization**
- **Multiscale and multitime elaboration and characterization technics.**
- **Modeling and simulation.**