



The metrological approach

A major key factor for the continuous improvement of the wood preservation laboratory

Luc MARTIN

Head , Quality & Metrology

Christine BAUDASSÉ, Nicolas LEMÉNAGER,
Marie-France THÉVENON



1

Presentation of the wood preservation laboratory



Research fields

- Rational use of lignocellulosic materials (wood)
- Service-life of ligno-cellulosic/wooden commodities
- Protection products & processes with low environmental impact

Main tasks

- Evaluation & comprehension of natural durability
- Development of products and wood processing with low environmental impact
- Evaluation of wood preservative chemicals



Activities

- Research
- Expertises
- Tests



Activities

- Research
- Expertises



- Tests → on durability of wood and wood-based products
→ on the protective efficacy of wood preservatives

Accreditation COFRAC for 8 tests

- **Termite tests** according to EN 117 and EN 118
- **Fungal tests** according to EN 113 and XP ENV 12038
Associated accelerated ageing procedures according to EN 73 and EN 84
- **Termite control products tests** according to XP X 41-550
Associated ageing procedure according to XP X 41-542





Determination of the resistance against wood-destroying agents !



Diversity of organisms



Evaluation of the impact
Biological degradation



Material heterogeneity



Determination of the resistance against wood-destroying agents !



Material heterogeneity

Evaluation of the impact
Biological degradation

Diversity of organisms




Staff Environment

Methods Material Equipments **Results**

Equipment pool

- 60 measurement devices including 14 measurement standards
- Management with SPLIT4® software
- Physical quantities: length, angle, weight, temperature, hygrometry, pressure, volume, air flow



2

Crucial steps
of EN 117 test
and
metrological
associated
issues



EN 117 standard

Wood preservatives - Determination of toxic values against *Reticulitermes* species (European termite)

Determination of the concentration at which the product completely prevents attack by termites of impregnated wood of a susceptible timber species

→ Scots Pine sapwood treated with a wood protection product using different concentrations



Method

- **Treatment** through a vacuum-pressure device (core impregnation)
- Measurement of the **quantity of product introduced** into the wood = product retention
- Drying (and eventually accelerated ageing)
- **Exposition of treated wood blocks** in culture flask to **250 termite workers** (8 weeks)
- Wood samples removed from the flasks
- Determination of the **number of survivors**
- **Visual rating of wood samples** (evaluation of the degree of attack)

Crucial steps

- Wood selection
- Counting the number of termites
- Visual quotation of the wood samples

Crucial steps

- **Wood selection**
- Counting the number of termites
- Visual quotation of the wood samples

Requirements

→ Reduce as much as possible the variability of the matrix

■ Timber to use

Exclusively sapwood containing little resin

■ Proportion of latewood < 30%

■ Ring orientation

Between 2.5 and 8 annual rings / cm


Contact angle > 10° to the broad faces of test specimens

■ Dimension of the blocks

Length: (50.0 ± 0.5) mm in longitudinal plan L

Width: (25.0 ± 0.5) mm in radial plan R



Thickness: (15.0 ± 0.5) mm tangential plan T



Metrological associated issue

■ Test gauge

Verification of orientation and proportion of annual rings




Help for operators

When the wood is cut: from tree log to plank and then to samples

When the blocks are selected for the test

After this first selection...

- Wood blocks must be measured



- Wood blocks must be weight

→ Determination of the density



Initial weight + density → Product retention calculation

Crucial steps

- Wood selection
- **Counting the number of termites**
- Visual quotation of the wood samples


Requirements

- Exposition to **250 termite workers** in a device
- For a complete evaluation: **40 devices → 10,000 termites!!**



Metrological associated issue

■ Termite-meter



Aspiration speed carefully chosen

→ Not to affect the survival rate or vitality of termites

Detection of the photo-sensible cell

→ Exact number given

→ Only termites are detected

Help for operators

Average time saving of **50%** for a test set up!!

Crucial steps

- Wood selection
- Counting the number of termites
- **Visual quotation of the wood samples**

Requirements

- After 8 weeks of biological exposure...



$(27 \pm 1)^{\circ} \text{C} > 75\% \text{HR}$

- ... Wood samples are removed
and **termites remained alive** are counting (by hand)
- To be declared valid, the **virulence controls** must present
 - a termite survival of 50%
 - a strong level of attack (level 4)

Expression of examination (1)

- Classification of attack by
 - its locality
 - its extent
 - its depth
 - Expression in accordance with the following schedule
 - Level 0: No attack
 - **Level 1: Attempted attack**
 - **Level 2: Slight attack**
 - Level 3: Average attack
 - Level 4: Strong attack
- Each level of attack can be a combinaison of 2 criteria

Expression of examination (2)

- Example → **Level 1** is reached for
 - (i) Superficial erosion of insufficient depth to be measured on an unlimited area of the test specimen
 - or (ii) Attack to a depth of 0.5 mm provided that this is restricted to an area or areas not more than 30 mm² in total
 - or (iii) A combination of (i) and (ii)



Metrological associated issue

1. Visual examination of each test specimen individually performed by 3 operators



2. Results commonly reviewed by all operators
3. When divergences occur, discussion leads to a given rating
→ « Laboratory result »

Inter-operator test according to ISO 5725-2

→ Operators qualification (not uncertainty estimation...)

3

Discussion
&
Conclusion

Zimbabwe

Reunion

Madagascar

Burkina Faso

Mexico

Laos

Brazil

EN 117 standard aims to know what is the **preventive protection level** against termites of a pine sapwood wood block containing a certain amount of product (**retention**) at a **given concentration**

1. For a **fixed concentration**, retention of wood preservative product in sample depends on **initial weight** and **volume**
 - **Wood selection** crucial step
 - **Test gauge** essential tool

Intermediary retention result estimated according to **ENV 13005**

2. Final result on level attack relies on **technical competence** of laboratory operators

Inter-operator comparison according to **ISO 5725-2**

Metrological point of view

- Methods such as **3D-image analysis** could substitute visual rating ?

Scientific point of view

- New way of measuring termite degradation should be developed according to **evolution of preservative products formulations** (toxicity + repellency)
 - Inter-laboratory cross-checking tests to develop new method based on **choice test**

Conclusion

- Accreditation is a guarantee of **laboratory's expertise** for research partners or to **respond to tenders** for European projects
- **Development of measuring instruments adapted to constraints of our testing has strengthened competence of operators and appears as an important factor in improving laboratory practices**
- Metrological approach is a **long-term process**, extension to
 - Other **standardized tests**
 - **Development of new methods**



cirad

AGRICULTURAL RESEARCH
FOR DEVELOPMENT

Working
together for
tomorrow's
agriculture

شكرا

**Thank you !
Merci !**

